

## Sommerhuber Wet Room Seat MONO

Wet room seat as module with 1° decline seat consisting of three Heat Storing Ceramics - backrest, seat and plinth.



Velvety soft ceramic surfaces with optimal hygiene aspects. Thanks to specific characteristics Heat Storing Ceramics emit mild areal infrared warmth, which enhances blood circulation and results in superior relaxation and maintaining of vitality.

Descriptive size: W/L/H ~72,6/ 62,9/ 133,9 cm or  
~ 72,6/ 48,9/ 135,5 cm

Article no.: 651 950

Delivery in separate components.  
Installation on site.

- Consisting of backrest, seat and plinth element.  
Decline seat by 1°
- standard width of joints ~ 3 mm
- recommended surface temperature: ~ 34°C +/- 2°C
- Optional: LED illumination



All measurements are nominal sizes.  
All measurements in cm.

This product is maintenance free.

For proper cleaning please use only cleaning agents with PH value between 5.0 - 9.0. Please use soft cleaning clothes only. Abrasive brushes or high-pressure cleaners could damage the glazed surface.